

HARVARD MEDICAL
SCHOOL



BRIGHAM AND WOMEN'S
HOSPITAL



Targeted Lipid Mediator Metabololipidomics LC-MS-MS

Eicosanoids and SPM workup

Center for Experimental Therapeutics and Reperfusion Injury,
Brigham and Women's Hospital, and
Harvard Medical School

Correspondence should be addressed to:

Prof. Charles N. Serhan
Hale Building for Transformative Medicine
60 Fenwood Road
Boston, MA 02115

Email: cnserhan@zeus.bwh.harvard.edu; Tel: 617-525-5001; Fax: 617-525-5017

We acknowledge the support of the NIH Grant P01GM095467

Lipid Mediator Metabolipidomics *LC-MS-MS workup*



© Serhan et al.

Commercially Available Standards

Cayman Catalog Numbers

Specialized Pro-Resolving Mediators

Lipoxin A ₄	(90410)
Lipoxin B ₄	(90420)
5S,15S-diHETE	(35280)
15epi-Lipoxin A ₄	(90415)
Resolvin D1	(10012554)
17epi-Resolvin D1	(13060)
Resolvin D2	(10007279)
Resolvin D3	(13834)
Resolvin D5	(10007280)
10(S),17(S)-DiHDHA (PDX)	(10008128)
Maresin 1	(10878)
Lipoxin A ₅	(10011453)

Biosynthesis Pathway Markers

AA	(90010)
5-HETE	(34210)
15(S)-HpETE	(44720)
5(S)-HpETE	(44230)
12-HETE	(34550)
15-HETE	(34700)
DHA	(90310)
17(S)-HpDHA	(13185)
14(S)-HDHA	(15253)
4-HDHA	(33200)
7-HDHA	(33300)
17-HDHA	(33650)
EPA	(90110)
5-HEPE	(32200)
12-HEPE	(32540)
15-HEPE	(32700)
18-HEPE	(32840)

AA bioactive Metabolome

PGE ₂	(14010)
PGD ₂	(12010)
PGF _{2α}	(16010)
20-OH-LTB ₄	(20190)
LTB ₄	(20110)
LTC ₄	(20210)
LTD ₄	(20310)
LTE ₄	(20410)
TXB ₂	(19030)





Table of Contents

Lipid Mediator Structures & Stereochemistry	4
Arachidonic Acid Metabolome	5
Eicosapentaenoic Acid Metabolome	10
Docosahexaenoic Acid Metabolome	14
Aspirin-Triggered Specialized Pro-resolving Mediators	19
Deuterium Labeled Internal Standards	23
Methods	25
Operationalization of lipid mediator metabololipidomics	27
High Pressure Liquid Chromatography	28
Mass Spectrometry	28
Methods validation	31
Linearity and attributed internal standard	31
Recovery, matrix interference and limit of detection	32
Repeatability and reproducibility	33
Chemical stability of selected endogenous lipid mediators identified in human serum	37
Multiple Reaction Monitoring Transitions	38
Lipid Mediators	39
Aspirin-Triggered Lipid Mediators	39
Further Metabolites of Lipid Mediators	40
Lipid Mediator Biosynthetic Pathways	41
Prostaglandins, Leukotrienes and Lipoxins	42
Leukotrienes	43
E-Series Resolvins	44
D-Series Resolvins	45
Aspirin-Triggered D-Series Resolvins	47
Protectins	48
Maresins	49

Lipid Mediator Metabolipidomics
LC-MS-MS workup





CNSerhanlab

CETRI

Center for Experimental Therapeutics
and Reperfusion Injury





Lipid Mediator Structures & Stereochemistry



CETRI



Center for Experimental Therapeutics
and Reperfusion Injury

Arachidonic Acid Metabolome

Abbreviation	Structure	Formula	Chemical Name	Trivial Name	Complete Stereochemistry
AA		$C_{20}H_{32}O_2$	Eicosatetraenoic acid	Arachidonic acid	5Z,8Z,11Z,14Z-Eicosatetraenoic acid
PGD₂		$C_{20}H_{32}O_5$	11-dehydro-prostaglandin F _{2α}	Prostaglandin D ₂	11-oxo-9α,15S-dihydroxy-5Z,13E-prostadienoic acid
PGE₂		$C_{20}H_{32}O_5$	9-dehydro-prostaglandin F _{2α}	Prostaglandin E ₂	9-oxo-11α,15S-dihydroxy-5Z,13E-prostadienoic acid
PGF_{2α}		$C_{20}H_{34}O_5$	9α,11α,15S-trihydroxy-prostaglandin F _{2α}	Prostaglandin F _{2α}	9α,11α,15S-trihydroxy-5Z,13E-prostadienoic acid

Abbreviation	Structure	Formula	Chemical Name	Trivial Name	Complete Stereochemistry
TxB₂		$C_{20}H_{34}O_6$	9 α ,11,15S-trihydroxy-thromboxane B ₂	Thromboxane B ₂	9 α ,11,15S-trihydroxy-5Z,13E-thrombadienoic acid
LTB₄		$C_{20}H_{32}O_4$	5S,12R-dihydroxy-LTB ₄	Leukotriene B ₄	5S,12R-dihydroxy-6Z,8E,10E,14Z-eicosatetraenoic acid
20-OH-LTB₄		$C_{20}H_{32}O_5$	5S,12R,20-trihydroxy-LTB ₄	20-hydroxy-Leukotriene B ₄	5S,12R,20-trihydroxy-6Z,8E,10E,14Z-eicosatetraenoic acid
20-COOH-LTB₄		$C_{20}H_{30}O_6$	5S,12R-dihydroxy-20-carboxy-LTB ₄	20-carboxy-Leukotriene B ₄	5S,12R-dihydroxy-6Z,8E,10E,14Z-20-carboxy-eicosatetraenoic acid


Abbreviation	Structure	Formula	Chemical Name	Trivial Name	Complete Stereochemistry
LXA₄		$C_{20}H_{32}O_5$	5S,6R,15S-trihydroxy-LXA ₄	Lipoxin A ₄	5S,6R,15S-trihydroxy-7E,9E,11Z,13E-eicosatetraenoic acid
LXB₄		$C_{20}H_{32}O_5$	5S,14R,15S-trihydroxy-LXB ₄	Lipoxin B ₄	5S,14R,15S-trihydroxy-6E,8Z,10E,12E-eicosatetraenoic acid
5S,15S-diHETE		$C_{20}H_{32}O_4$	5S,15S-dihydroxy-ETE	5S,15S-dihydroxy-eicosatetraenoic acid	5S,15S-dihydroxy-6E,8Z,11Z,13E-eicosatetraenoic acid
5-HETE		$C_{20}H_{32}O_3$	5-hydroxy-ETE	5-hydroxy-eicosatetraenoic acid	5-hydroxy-6E,8Z,11Z,14Z-eicosatetraenoic acid

Abbreviation	Structure	Formula	Chemical Name	Trivial Name	Complete Stereochemistry
12-HETE		$C_{20}H_{32}O_3$	12-hydroxy-ETE	12-hydroxy-eicosatetraenoic acid	12-hydroxy-5Z,8Z,10E,14Z-eicosatetraenoic acid
15-HETE		$C_{20}H_{32}O_3$	15-hydroxy-ETE	15-hydroxy-eicosatetraenoic acid	15-hydroxy-5Z,8Z,11Z,13E-eicosatetraenoic acid

Lipid Mediator Metabololipidomics ***LC-MS-MS workup***



© Serhan et al.











CNSerhanlab




CETRI



Center for Experimental Therapeutics
and Reperfusion Injury

Eicosapentaenoic Acid Metabolome

Abbreviation	Structure	Formula	Chemical Name	Trivial Name	Complete Stereochemistry
EPA		$C_{20}H_{30}O_2$	Eicosapentaenoic acid	Eicosapentaenoic acid	5Z,8Z,11Z,14Z,17Z-eicosapentaenoic acid
LXA₅		$C_{20}H_{30}O_5$	5S,6R,15S-trihydroxy-LXA ₅	Lipoxin A ₅	5S,6R,15S-trihydroxy-7E,9E,11Z,13E,17Z-eicosapentaenoic acid
LXB₅		$C_{20}H_{30}O_5$	5S,14R,15S-trihydroxy-LXB ₅	Lipoxin B ₅	5S,14R,15S-trihydroxy-6E,8Z,10E,12E,17Z-eicosapentaenoic acid
5S,15S-diHEPE		$C_{20}H_{30}O_4$	5S,15S-dihydroxy-EPE	5S,15S-dihydroxy-eicosapentaenoic acid	5S,15S-dihydroxy-6E,8Z,11Z,13E,17Z-eicosapentaenoic acid

Abbreviation	Structure	Formula	Chemical Name	Trivial Name	Complete Stereochemistry
RVE1		$C_{20}H_{30}O_5$	5 <i>S</i> ,12 <i>R</i> ,18 <i>R</i> -trihydroxy-EPE	Resolvin E1	5 <i>S</i> ,12 <i>R</i> ,18 <i>R</i> -trihydroxy-6 <i>Z</i> ,8 <i>E</i> ,10 <i>E</i> ,14 <i>Z</i> ,16 <i>E</i> -eicosapentaenoic acid
RVE2		$C_{20}H_{30}O_4$	5 <i>S</i> ,18 <i>R</i> -dihydroxy-EPE	Resolvin E2	5 <i>S</i> ,18 <i>R</i> -dihydroxy-6 <i>E</i> ,8 <i>Z</i> ,11 <i>Z</i> ,14 <i>Z</i> ,16 <i>E</i> -eicosapentaenoic acid
RVE3		$C_{20}H_{30}O_4$	17 <i>R</i> ,18 <i>R</i> -dihydroxy-EPE	Resolvin E3	17 <i>R</i> ,18 <i>R</i> -dihydroxy-5 <i>Z</i> ,8 <i>Z</i> ,11 <i>Z</i> ,13 <i>E</i> ,15 <i>E</i> -eicosapentaenoic acid
5-HEPE		$C_{20}H_{30}O_3$	5-hydroxy-EPE	5-hydroxy-eicosapentaenoic acid	5-hydroxy-6 <i>E</i> ,8 <i>Z</i> ,11 <i>Z</i> ,14 <i>Z</i> ,17 <i>Z</i> -eicosapentaenoic acid

Abbreviation	Structure	Formula	Chemical Name	Trivial Name	Complete Stereochemistry
12-HEPE		$C_{20}H_{30}O_3$	12-hydroxy-EPE	12-hydroxy-eicosapentaenoic acid	12-hydroxy-5Z,8Z,10E,14Z,17Z-eicosapentaenoic acid
15-HEPE		$C_{20}H_{30}O_3$	15-hydroxy-EPE	15-hydroxy-eicosapentaenoic acid	15-hydroxy-5Z,8Z,11Z,13E,17Z-eicosapentaenoic acid
18-HEPE		$C_{20}H_{30}O_3$	18-hydroxy-EPE	18-hydroxy-eicosapentaenoic acid	18-hydroxy-5Z,8Z,11Z,14Z,16E-eicosapentaenoic acid











CNSerhanlab

CETRI



Center for Experimental Therapeutics
and Reperfusion Injury

Docosaehaenoic Acid Metabolome


Abbreviation	Structure	Formula	Chemical Name	Trivial Name	Complete Stereochemistry
DHA		$C_{22}H_{32}O_2$	Docosahexaenoic acid	Docosahexaenoic acid	4Z,7Z,10Z,13Z,16Z,19Z-docosahexaenoic acid
RvD1		$C_{22}H_{32}O_5$	7S,8R,17S-trihydroxy-DHA	Resolvin D1	7S,8R,17S-trihydroxy-4Z,9E,11E,13Z,15E,19Z-docosahexaenoic acid
RvD2		$C_{22}H_{32}O_5$	7S,16R,17S-trihydroxy-DHA	Resolvin D2	7S,16R,17S-trihydroxy-4Z,8E,10Z,12E,14E,19Z-docosahexaenoic acid
RvD3		$C_{22}H_{32}O_5$	4S,11R,17S-trihydroxy-DHA	Resolvin D3	4S,11R,17S-trihydroxy-5Z,7E,9E,13Z,15E,19Z-docosahexaenoic acid

Abbreviation	Structure	Formula	Chemical Name	Trivial Name	Complete Stereochemistry
RvD4		$C_{22}H_{32}O_5$	4S,5R,17S-trihydroxy-DHA	Resolvin D4	4S,5R,17S-trihydroxy-docosahexaenoic acid
RvD5		$C_{22}H_{32}O_4$	7S,17S-dihydroxy-DHA	Resolvin D5	7S,17S-dihydroxy-docosahexaenoic acid
RvD6		$C_{22}H_{32}O_4$	4S,17S-dihydroxy-DHA	Resolvin D6	4S,17S-dihydroxy-docosahexaenoic acid
MaR1		$C_{22}H_{32}O_4$	7R,14S-dihydroxy-DHA	Maresin 1	7R,14S-dihydroxy-docosahexaenoic acid

Abbreviation	Structure	Formula	Chemical Name	Trivial Name	Complete Stereochemistry
4S,14S-diHDHA		$C_{22}H_{32}O_4$	4S,14S-dihydroxy-DHA	4,14-diHDHA	4S,14S-dihydroxydocosahexaenoic acid
PD1		$C_{22}H_{32}O_4$	10R,17S-dihydroxy-DHA	[Neuro] Protectin D1	10R,17S-dihydroxydocosahexaenoic acid
4-HDHA		$C_{22}H_{32}O_3$	4-hydroxy-DHA	4-hydroxy-DHA	4-hydroxydocosahexaenoic acid
7-HDHA		$C_{22}H_{32}O_3$	7-hydroxy-DHA	7-hydroxy-DHA	7-hydroxydocosahexaenoic acid

Abbreviation	Structure	Formula	Chemical Name	Trivial Name	Complete Stereochemistry
14-HDHA		$C_{22}H_{32}O_3$	14-hydroxy-DHA	14-hydroxy-DHA	14-hydroxy-4Z,7Z,10Z,12E,16Z,19Z-docosahexaenoic acid
17-HDHA		$C_{22}H_{32}O_3$	17-hydroxy-DHA	17-hydroxy-DHA	17-hydroxy-4Z,7Z,10Z,13Z,15E,19Z-docosahexaenoic acid





Lipid Mediator Metabolipidomics
LC-MS-MS workup








© Serhan et al.




Aspirin-Triggered Specialized Pro-resolving Mediators

Abbreviation	Structure	Formula	Chemical Name	Trivial Name	Complete Stereochemistry
AT-LXA₄		$C_{20}H_{32}O_5$	5S,6R,15R-trihydroxy-LXA ₄	15-epi-Lipoxin A ₄	5S,6R,15R-trihydroxy-7E,9E,11Z,13E-eicosatetraenoic acid
AT-LXB₄		$C_{20}H_{32}O_5$	5S,14R,15R-trihydroxy-LXB ₄	15-epi-Lipoxin B ₄	5S,14R,15R-trihydroxy-6E,8Z,10E,12E-eicosatetraenoic acid
AT-RvD1		$C_{22}H_{32}O_5$	7S,8R,17R-trihydroxy-DHA	Aspirin-triggered Resolvin D1	7S,8R,17R-trihydroxy-4Z,9E,11E,13Z,15E,19Z-docosahexaenoic acid
AT-RvD2		$C_{22}H_{32}O_5$	7S,16R,17R-trihydroxy-DHA	Aspirin-triggered Resolvin D2	7S,16R,17R-trihydroxy-4Z,8E,10Z,12E,14E,19Z-docosahexaenoic acid

Abbreviation	Structure	Formula	Chemical Name	Trivial Name	Complete Stereochemistry
AT-RvD3		$C_{22}H_{32}O_5$	4 <i>S</i> ,11 <i>R</i> ,17 <i>R</i> -trihydroxy-DHA	Aspirin-triggered Resolvin D3	4 <i>S</i> ,11 <i>R</i> ,17 <i>R</i> -trihydroxy-5 <i>Z</i> ,7 <i>E</i> ,9 <i>E</i> ,13 <i>Z</i> ,15 <i>E</i> ,19 <i>Z</i> -docosahexaenoic acid
AT-RvD4		$C_{22}H_{32}O_5$	4 <i>S</i> ,5,17 <i>R</i> -trihydroxy-DHA	Aspirin-triggered Resolvin D4	4 <i>S</i> ,5 <i>R</i> ,17 <i>R</i> -trihydroxy-6 <i>E</i> ,8 <i>E</i> ,10 <i>Z</i> ,13 <i>Z</i> ,15 <i>E</i> ,19 <i>Z</i> -docosahexaenoic acid
AT-RvD5		$C_{22}H_{32}O_4$	7 <i>S</i> ,17 <i>R</i> -dihydroxy-DHA	Aspirin-triggered Resolvin D5	7 <i>S</i> ,17 <i>R</i> -dihydroxy-4 <i>Z</i> ,8 <i>E</i> ,10 <i>Z</i> ,13 <i>Z</i> ,15 <i>E</i> ,19 <i>Z</i> -docosahexaenoic acid
AT-RvD6		$C_{22}H_{32}O_4$	4 <i>S</i> ,17 <i>R</i> -dihydroxy-DHA	Aspirin-triggered Resolvin D6	4 <i>S</i> ,17 <i>R</i> -dihydroxy-5 <i>E</i> ,7 <i>Z</i> ,10 <i>Z</i> ,13 <i>Z</i> ,15 <i>E</i> ,19 <i>Z</i> -docosahexaenoic acid

Abbreviation	Structure	Formula	Chemical Name	Trivial Name	Complete Stereochemistry
AT-PD1		$C_{22}H_{32}O_4$	10 <i>R</i> ,17 <i>R</i> -dihydroxy-DHA	Aspirin-triggered [Neuro] Protectin D1	10 <i>R</i> ,17 <i>R</i> -dihydroxy-4 <i>Z</i> ,7 <i>Z</i> ,11 <i>E</i> ,13 <i>E</i> ,15 <i>Z</i> ,19 <i>Z</i> -docosahexaenoic acid





Lipid Mediator Metabolipidomics
LC-MS-MS workup




CETRI


Center for Experimental Therapeutics
and Reperfusion Injury

Deuterium Labeled Internal Standards


Abbreviation	Structure	Formula	Chemical Name	Trivial Name	Complete Stereochemistry
d₄-PGE₂		$C_{20}H_{28}D_4O_5$	9-dehydro- prostaglandin F _{2α} -d ₄	d ₄ -Prostaglandin E ₂	9-oxo- 11α,15S- dihydroxy- 5Z,13E- 3,3,4,4-d ₄ - prostaglandienoic acid
d₅-LXA₄		$C_{20}H_{27}D_5O_5$	5S,6R,15S- trihydroxy- LXA ₄ -d ₅	d ₅ -Lipoxin A ₄	5S,6R,15S- trihydroxy- 7E,9E,11Z,13E- 19,19,20,20,20-d ₅ - eicosatetraenoic acid
d₄-LTB₄		$C_{20}H_{28}D_4O_4$	5S,12R- dihydroxy- LTB ₄ -d ₄	d ₄ -Leukotriene B ₄	5S,12R- dihydroxy- 6Z,8E,10E,14Z- 6,7,14,15-d ₄ - eicosatetraenoic acid
d₈-5S-HETE		$C_{20}H_{24}D_8O_3$	5S- hydroxy- ETE-d ₈	d ₈ -5S- hydroxy- eicosatetraenoic acid	5S- hydroxy- 6E,8Z,11Z,14Z- 5,6,8,9,11,12,14,15-d ₈ - eicosatetraenoic acid

Abbreviation	Structure	Formula	Chemical Name	Trivial Name	Complete Stereochemistry
d₅-RvD2		$C_{22}H_{27}D_5O_5$	7 <i>S</i> ,16 <i>R</i> ,17 <i>S</i> -trihydroxy-DHA -d ₅	d ₅ -Resolvin D2	7 <i>S</i> ,16 <i>R</i> ,17 <i>S</i> -trihydroxy-4 <i>Z</i> ,8 <i>E</i> ,10 <i>Z</i> ,12 <i>E</i> ,14 <i>E</i> ,19 <i>Z</i> -19,19,20,20,20,20-d ₅ -docosahexaenoic acid

Lipid Mediator Metabololipidomics ***LC-MS-MS workup***



© Serhan et al.



CNSerhanlab

CETRI


Center for Experimental Therapeutics
and Reperfusion Injury

Methods

Lipid Mediator Metabololipidomics LC-MS-MS workup



Operationalization of lipid mediator metabololipidomics



Identification and signature profiles for pro-resolving and inflammatory lipid mediators in human tissue
 Colas RA, Shinohara M, Dalli J, Chiang N, Serhan CN. Am J Physiol Cell Physiol. 2014

Lipid Mediator Metabolipidomics

LC-MS-MS workup



© Serhan et al.

Shimadzu LC20AD HPLC Settings		
Solvent A	H ₂ O (0.01% acetic acid)	
Solvent B	Methanol (0.01% acetic acid)	
Flow rate	0.4mL/min	
Column	Eclipse plus C18, 1.8 μm, 4.6 x 100 mm, AGILENT	
Column Heater	ThermaSphere TS-130, 50°C, PHENOMENEX	
Gradient	Time (min)	Solvent (%B)
	0.0	55
	2.0	55
	10.0	85
	18.0	98
	20.0	98
	20.1	55
	24.5	End

Table 1: High Performance Liquid Chromatography settings.

AB Sciex 5500 Q TRAP MS setting for MRM and EPI		
	MRM	EPI
Curtain Gas	25	25
Collision Gas	Medium	Medium
Ion Spray Voltage	-4000	-4000
Temperature	500 °C	500 °C
Ion Source Gas 1	40	40
Ion Source Gas 2	40	40

Table 2: Mass Spectrometer settings for Multiple Reaction Monitoring and Enhance Product Ion for AB Sciex 5500 Q TRAP.

AB Sciex 6500 Q TRAP MS setting for MRM and EPI		
	MRM	EPI
Curtain Gas	30	30
Collision Gas	Medium	Medium
Ion Spray Voltage	-4500	-4500
Temperature	580 °C	580 °C
Ion Source Gas 1	80	80
Ion Source Gas 2	70	70

Table 3: Mass Spectrometer settings for Multiple Reaction Monitoring and Enhance Product Ion for AB Sciex 6500 Q TRAP.

Lipid Mediator Metabololipidomics LC-MS-MS workup



© Serhan et al.

AB Sciex 5500 Q TRAP				
Individual lipid mediator, pathway marker, and further metabolite setting for MRM (V)				
Compound	DP	EP	CE	CXP
AA	-90	-10	-41	-22
PGD ₂	-80	-10	-30	-15
PGE ₂	-80	-10	-30	-15
PGF _{2α}	-80	-10	-30	-15
TxB ₂	-80	-10	-30	-15
LTB ₄	-110	-10	-22	-17
20-OH-LTB ₄	-100	-10	-26	-17
20-COOH-LTB ₄	-90	-10	-26	-17
LXA ₄	-95	-10	-22	-11
LXB ₄	-90	-10	-30	-21
5,15-diHETE	-110	-10	-22	-17
5-HETE	-90	-10	-21	-22
12-HETE	-70	-10	-21	-15
15-HETE	-65	-10	-19	-18
EPA	-90	-10	-41	-22
LXA ₅	-95	-10	-22	-11
LXB ₅	-90	-10	-30	-21
5,15-diHEPE	-110	-10	-22	-17
RvE1	-80	-10	-24	-16
RvE2	-40	-10	-22	-17
RvE3	-40	-10	-22	-17
5-HEPE	-60	-10	-22	-11
12-HEPE	-70	-10	-21	-15
15-HEPE	-65	-10	-19	-8
18-HEPE	-65	-10	-18	-23
DHA	-90	-10	-41	-22
RvD1	-90	-10	-22	-13
RvD2	-90	-10	-24	-13
RvD3	-80	-10	-28	-19
RvD5	-40	-10	-22	-17
RvD6	-40	-10	-22	-17
MaR1	-40	-10	-22	-17
4,14-diHDHA	-40	-10	-22	-17
PD1	-40	-10	-22	-17
4-HDHA	-40	-10	-18	-13
7-HDHA	-40	-10	-18	-13
14-HDHA	-70	-10	-21	-15
17-HDHA	-80	-10	-20	-18
d ₄ -PGE ₂	-80	-10	-30	-11
d ₄ -LXA ₄	-135	-10	-24	-17
d ₄ -LTB ₄	-135	-10	-24	-17
d ₈ -5S-HETE	-100	-10	-22	-11
d ₅ -RvD2	-90	-10	-24	-13

Table 4: Mass Spectrometer settings for Multiple Reaction Monitoring for AB Sciex 5500 Q TRAP. DP=Declustering Potential, EP=Entrance potential, CE=Collision Energy, CXP=Collision Cell Exit Potential

Lipid Mediator Metabololipidomics LC-MS-MS workup



© Serhan et al.

AB Sciex 6500 Q TRAP				
Individual lipid mediator, pathway marker, and further metabolite setting for MRM (V)				
Compound	DP	EP	CE	CXP
AA	-10.0	-10.0	-10.0	-12.0
PGD ₂	-30.0	-10.0	-22.0	-8.0
PGE ₂	-30.0	-10.0	-22.0	-8.0
PGF _{2α}	-60.0	-10.0	-33.3	-12.0
TxB ₂	-70.0	-10.0	-25.0	-22.0
LTB ₄	-45.0	-10.0	-21.0	-10.8
20-OH-LTB ₄	-60.0	-10.0	-22.0	-8.0
20-COOH-LTB ₄	-60.0	-10.0	-22.0	-8.0
LXA ₄	-30.0	-10.0	-18.2	-10.0
LXB ₄	-45.0	-10.0	-21.5	-12.0
5,15-diHETE	-35.0	-10.0	-18.0	-12.0
5-HETE	-30.0	-10.0	-18.0	-9.5
12-HETE	-30.0	-10.0	-18.0	-9.5
15-HETE	-30.0	-10.0	-18.0	-9.5
EPA	-10.0	-10.0	-10.0	-12.0
LXA ₅	-30.0	-10.0	-18.6	-10.0
LXB ₅	-30.0	-10.0	-21.2	-10.0
5,15-diHEPE	-30.0	-12.0	-17.0	-12.0
RvE1	-50.0	-10.0	-21.1	-10.0
RvE2	-30.0	-12.0	-18.6	-10.0
RvE3	-35.0	-10.0	-20.5	-14.0
5-HEPE	-30.0	-10.0	-18.0	-9.5
12-HEPE	-30.0	-10.0	-18.0	-9.5
15-HEPE	-30.0	-10.0	-18.0	-9.5
18-HEPE	-30.0	-10.0	-18.0	-9.5
DHA	-10.0	-10.0	-10.0	-12.0
RvD1	-40.0	-10.0	-19.0	-9.0
RvD2	-40.0	-10.0	-27.8	-10.0
RvD3	-55.0	-10.0	-24.6	-9.5
RvD5	-42.0	-10.0	-18.0	-11.0
RvD6	-50.0	-10.0	-20.0	-11.0
MaR1	-38.0	-10.0	-19.6	-11.0
4,14-diHDHA	-50.0	-10.0	-20.0	-11.0
PD1	-45.0	-10.0	-20.8	-9.5
4-HDHA	-30.0	-10.0	-18.0	-9.5
7-HDHA	-30.0	-10.0	-18.0	-9.5
14-HDHA	-30.0	-10.0	-18.0	-9.5
17-HDHA	-30.0	-10.0	-18.0	-9.5
d ₄ -PGE ₂	-30.0	-10.0	-26.0	-14.0
d ₄ -LXA ₄	-40.0	-10.0	-17.5	-7.5
d ₄ -LTB ₄	-55.0	-10.0	-21.7	-10.0
d ₈ -5S-HETE	-40.0	-10.0	-24.0	-12.0
d ₅ -RvD2	-40.0	-10.0	-20.3	-13.0

Table 5: Mass Spectrometer settings for Multiple Reaction Monitoring for AB Sciex 5500 Q TRAP. DP=Declustering Potential, EP=Entrance potential, CE=Collision Energy, CXP=Collision Cell Exit Potential

Lipid Mediator Metabolipidomics LC-MS-MS workup




© Serhan et al.

Compound name	Linear equation	r ²	Internal standard
RvD1	y=0.00002081x	0.999	d ₅ -RvD2
RvD2	y=0.00009118x	0.997	d ₅ -RvD2
RvD3	y=0.000008163x	0.990	d ₅ -RvD2
RvD5	y=0.00001137x	0.993	d ₄ -LTB ₄
PD1	y=0.00001425x	0.995	d ₄ -LTB ₄
MaR1	y=0.00007534x	0.999	d ₄ -LTB ₄
17-HDHA	y=0.00001132x	0.996	d ₈ -5HETE
14-HDHA	y=0.00001916x	0.999	d ₈ -5HETE
7-HDHA	y=0.00001096x	0.996	d ₈ -5HETE
4-HDHA	y=0.00001040x	0.996	d ₈ -5HETE
LXA ₅	y=0.00002043x	0.999	d ₅ -LXA ₄
RvE1	y=0.00004358x	0.999	d ₅ -LXA ₄
RvE2	y=0.0001647x	0.996	d ₄ -LTB ₄
RvE3	y=0.0001047x	0.998	d ₄ -LTB ₄
18-HEPE	y=0.00002131x	0.998	d ₈ -5HETE
15-HEPE	y=0.00002002x	0.999	d ₈ -5HETE
12-HEPE	y=0.00001899x	0.997	d ₈ -5HETE
5-HEPE	y=0.00001301x	0.998	d ₈ -5HETE
LXA ₄	y=0.00001607x	0.996	d ₅ -LXA ₄
LXB ₄	y=0.00007325x	0.995	d ₅ -LXA ₄
5S,15S-diHETE	y=0.00008011x	0.997	d ₄ -LTB ₄
LTB ₄	y=0.00001074x	0.991	d ₄ -LTB ₄
20-OH-LTB ₄	y=0.00001524x	0.999	d ₄ -LTB ₄
20-COOH-LTB ₄	y=0.00007931x	0.999	d ₄ -LTB ₄
PGD ₂	y=0.00001019x	0.993	d ₄ -PGE ₂
PGE ₂	y=0.00001698x	0.998	d ₄ -PGE ₂
PGF _{2α}	y=0.00001643x	0.997	d ₄ -PGE ₂
TXB ₂	y=0.00001054x	0.997	d ₄ -PGE ₂
15-HETE	y=0.00001377x	0.999	d ₈ -5HETE
12-HETE	y=0.00001142x	0.994	d ₈ -5HETE
5-HETE	y=0.00001010x	0.990	d ₈ -5HETE

Table 6: Linear equation, coefficient of determination (r²) and attributed internal standard

Lipid Mediator Metabololipidomics

LC-MS-MS workup



© Serhan et al.

	Recovery (%) in plasma	Interference (% signal remaining)	LOD 5500 (pg) in		LOD 6500 (pg) in	
			phase	plasma	phase	plasma
Internal standards						
d ₅ -RvD2	73 ± 16	28 ± 9	0.17	0.43	0.04	0.09
d ₅ -LXA ₄	75 ± 2	7 ± 2	0.10	0.25	0.02	0.05
d ₄ -LTB ₄	80 ± 3	18 ± 2	0.10	0.51	0.05	0.22
d ₄ -PGE ₂	63 ± 13	21 ± 6	0.19	0.47	0.09	0.23
d ₈ -5S-HETE	68 ± 3	28 ± 5	0.10	0.26	0.07	0.15
d₅-authentic lipid mediators						
d ₅ -10S,17S-diHDHA	88 ± 4	8 ± 2	0.06	0.45	0.03	0.15
d ₅ -5S,15S-diHEPE	94 ± 6	9 ± 5	0.13	0.95	0.03	0.15
d ₅ -15-HEPE	63 ± 12	26 ± 4	0.09	0.23	0.10	0.19

Table 7: Determination of recovery, interference and LOD. Plasma recoveries, interference and limit of detection (LOD) of the internal standards and authentic d₅-lipid mediators used as representative of specific chromatographic regions and physical properties of compounds. LOD was assessed in phase and in the matrix. Mean ± SEM, n = 6 for recoveries and interference, d = 3 for LOD.

Lipid Mediator Metabololipidomics LC-MS-MS workup




© Serhan et al.

Compound name	Nominal concentration (pg)	Intraday (repeatability)				Interday (reproducibility)			
		Measured (pg)	%RSD	Accuracy	Measured (pg)	%RSD	Accuracy		
RvD1	5	5.0 ± 0.1	1%	99%	4.6 ± 0.1	2%	92%		
	20	17.2 ± 0.7	7%	86%	23.2 ± 0.7	5%	84%		
	100	109.8 ± 2.2	4%	90%	118.5 ± 1.8	3%	82%		
	200	227.8 ± 2.5	2%	86%	233.5 ± 1.8	2%	83%		
RvD2	5	5.2 ± 0.1	4%	96%	4.6 ± 0.1	4%	93%		
	20	16.9 ± 0.7	7%	85%	22.6 ± 0.5	4%	87%		
	100	113.0 ± 2.3	4%	87%	116.5 ± 1.6	3%	84%		
	200	227.8 ± 3.1	3%	86%	235.5 ± 2.5	2%	82%		
RvD3	5	5.7 ± 0.1	3%	86%	5.4 ± 0.1	4%	91%		
	20	19.8 ± 0.5	4%	99%	26.0 ± 0.3	2%	70%		
	100	118.3 ± 1.4	2%	82%	125.3 ± 0.9	2%	75%		
	200	226.8 ± 2.2	2%	87%	226.8 ± 1.7	1%	87%		
RvD5	5	4.8 ± 0.1	5%	97%	4.6 ± 0.04	2%	93%		
	20	16.7 ± 0.6	6%	83%	24.0 ± 0.3	2%	80%		
	100	114.3 ± 0.8	1%	86%	117.0 ± 0.0	0%	83%		
	200	217.0 ± 2.6	2%	92%	217.0 ± 2.4	2%	92%		
PD1	5	4.6 ± 0.1	5%	93%	4.4 ± 0.05	2%	88%		
	20	16.7 ± 0.3	3%	84%	23.7 ± 0.5	4%	82%		
	100	114.0 ± 0.7	1%	86%	118.5 ± 1.0	2%	82%		
	200	218.5 ± 1.3	1%	91%	220.0 ± 1.3	1%	90%		
MaR1	5	4.6 ± 0.2	9%	91%	4.2 ± 0.02	1%	85%		
	20	16.1 ± 0.1	1%	81%	22.3 ± 0.5	4%	89%		
	100	107.0 ± 3.5	7%	93%	110.8 ± 1.4	3%	89%		
	200	221.3 ± 2.3	2%	89%	229.0 ± 1.8	2%	86%		
17-HDHA	5	3.5 ± 0.1	4%	69%	4.4 ± 0.05	2%	87%		
	20	17.8 ± 0.3	3%	89%	23.7 ± 0.8	6%	82%		
	100	115.0 ± 1.1	2%	85%	120.0 ± 0.6	1%	80%		
	200	222.0 ± 1.2	1%	89%	227.5 ± 1.6	1%	86%		
14-HDHA	5	3.2 ± 0.1	5%	64%	4.3 ± 0.02	1%	86%		
	20	17.3 ± 0.4	4%	86%	22.2 ± 0.4	3%	89%		
	100	114.8 ± 2.4	4%	85%	109.5 ± 2.3	4%	91%		
	200	212.3 ± 2.3	2%	94%	233.5 ± 0.9	1%	83%		
7-HDHA	5	3.5 ± 0.1	4%	70%	4.5 ± 0.1	3%	90%		
	20	19.7 ± 0.1	1%	99%	25.8 ± 0.4	3%	71%		
	100	123.8 ± 2.3	4%	76%	128.0 ± 2.5	4%	72%		
	200	236.5 ± 1.0	1%	82%	240.5 ± 2.4	2%	80%		
4-HDHA	5	3.6 ± 0.1	6%	73%	4.6 ± 0.02	1%	91%		
	20	19.8 ± 0.3	2%	99%	26.0 ± 0.3	2%	70%		
	100	125.8 ± 2.0	3%	74%	128.5 ± 1.6	2%	72%		
	200	236.3 ± 2.1	2%	82%	240.3 ± 1.4	1%	80%		

Table 8 (part 1/4): Precision and accuracy of LM metabololipidomics. %RSD, relative standard deviation; %RSD = (SEM/mean) x 100; mean ± SEM; n = 4

Lipid Mediator Metabololipidomics LC-MS-MS workup




© Serhan et al.

Compound name	Nominal concentration (pg)	Intraday (repeatability)			Interday (reproducibility)		
		Measured (pg)	%RSD	Accuracy	Measured (pg)	%RSD	Accuracy
LXA ₅	5	5.6 ± 0.1	2%	87%	5.0 ± 0.1	4%	99%
	20	18.0 ± 0.3	3%	90%	21.6 ± 0.7	6%	92%
	100	111.0 ± 2.6	5%	89%	110.3 ± 1.9	3%	90%
	200	212.3 ± 2.8	3%	94%	240.3 ± 3.3	3%	80%
RvE1	5	5.8 ± 0.1	4%	85%	5.1 ± 0.1	3%	99%
	20	17.4 ± 0.1	1%	87%	24.3 ± 0.3	2%	78%
	100	113.3 ± 2.3	4%	87%	121.8 ± 1.7	3%	78%
	200	242.0 ± 3.7	3%	79%	241.8 ± 1.7	1%	79%
RvE2	5	5.5 ± 0.4	14%	91%	4.7 ± 0.1	3%	93%
	20	16.4 ± 0.7	8%	82%	20.7 ± 0.7	6%	97%
	100	108.8 ± 2.8	5%	91%	101.2 ± 2.9	6%	99%
	200	199.3 ± 3.3	3%	100%	240.0 ± 1.1	1%	80%
RvE3	5	6.1 ± 0.3	8%	78%	5.0 ± 0.4	16%	99%
	20	17.8 ± 1.7	17%	89%	23.7 ± 1.6	12%	81%
	100	105.0 ± 3.0	6%	95%	116.5 ± 1.0	2%	84%
	200	218.5 ± 3.4	3%	91%	227.8 ± 3.1	3%	86%
18-HEPE	5	4.1 ± 0.1	3%	82%	4.3 ± 0.03	2%	87%
	20	17.6 ± 0.2	2%	88%	23.3 ± 0.5	4%	83%
	100	117.5 ± 2.0	3%	83%	121.3 ± 1.5	2%	79%
	200	239.0 ± 2.3	2%	81%	241.8 ± 1.5	1%	79%
15-HEPE	5	4.3 ± 0.1	3%	87%	4.4 ± 0.04	2%	88%
	20	17.6 ± 0.2	2%	88%	22.8 ± 0.5	4%	86%
	100	115.5 ± 2.5	4%	85%	122.3 ± 1.1	2%	78%
	200	231.3 ± 1.7	1%	84%	239.3 ± 1.7	1%	80%
12-HEPE	5	4.6 ± 0.04	2%	93%	4.4 ± 0.1	3%	88%
	20	18.1 ± 0.2	1%	91%	22.6 ± 0.4	3%	87%
	100	122.5 ± 1.3	2%	78%	117.3 ± 1.7	3%	83%
	200	231.3 ± 2.3	2%	84%	227.8 ± 1.3	1%	86%
5-HEPE	5	4.6 ± 0.05	2%	91%	4.6 ± 0.04	2%	92%
	20	18.1 ± 0.2	2%	91%	23.4 ± 0.5	4%	83%
	100	115.3 ± 1.3	2%	85%	119.5 ± 0.5	1%	81%
	200	223.5 ± 2.0	2%	88%	233.3 ± 1.0	1%	83%

Table 8 (part 2/4): Precision and accuracy of LM metabololipidomics. %RSD, relative standard deviation; %RSD = (SEM/mean) x 100; mean ± SEM; n = 4

Lipid Mediator Metabololipidomics LC-MS-MS workup



© Serhan et al.

Compound name	Nominal concentration (pg)	Intraday (repeatability)			Interday (reproducibility)		
		Measured (pg)	%RSD	Accuracy	Measured (pg)	%RSD	Accuracy
LXA ₄	5	5.5 ± 0.1	2%	90%	4.9 ± 0.1	4%	98%
	20	17.7 ± 0.8	8%	88%	23.5 ± 0.1	1%	82%
	100	114.0 ± 2.1	4%	86%	119.0 ± 1.1	2%	81%
	200	226.3 ± 2.1	2%	87%	233.5 ± 1.5	1%	83%
LXB ₄	5	6.2 ± 0.1	4%	76%	5.4 ± 0.1	3%	92%
	20	18.4 ± 0.6	6%	92%	22.9 ± 0.4	3%	85%
	100	116.8 ± 1.7	3%	83%	114.5 ± 1.2	2%	86%
	200	220.0 ± 4.5	4%	90%	231.8 ± 1.6	1%	84%
5S,15S-diHETE	5	5.1 ± 0.1	4%	98%	5.4 ± 0.1	2%	93%
	20	20.3 ± 0.5	4%	98%	26.4 ± 0.2	1%	68%
	100	123.3 ± 1.4	2%	77%	129.0 ± 2.0	3%	71%
	200	246.5 ± 3.7	3%	77%	251.0 ± 6.1	5%	75%
LTB ₄	5	5.0 ± 0.03	1%	99%	5.8 ± 0.1	2%	84%
	20	22.3 ± 0.3	2%	89%	27.5 ± 0.4	3%	63%
	100	125.0 ± 0.6	1%	75%	130.0 ± 1.8	3%	70%
	200	231.0 ± 3.2	3%	85%	226.8 ± 2.9	3%	87%
20-OH-LTB ₄	5	5.4 ± 0.1	4%	93%	4.9 ± 0.03	1%	98%
	20	17.1 ± 0.5	6%	86%	22.6 ± 0.4	3%	87%
	100	111.3 ± 2.7	5%	89%	116.0 ± 0.8	1%	84%
	200	223.3 ± 1.1	1%	88%	229.0 ± 3.0	3%	86%
20-COOH-LTB ₄	5	5.1 ± 0.4	16%	99%	5.2 ± 0.1	2%	97%
	20	18.0 ± 0.9	9%	90%	21.4 ± 1.6	13%	93%
	100	117.5 ± 1.3	2%	83%	114.0 ± 2.2	4%	86%
	200	227.3 ± 2.3	2%	86%	253.8 ± 3.4	3%	73%
PGD ₂	5	5.7 ± 0.1	2%	87%	5.2 ± 0.1	4%	96%
	20	17.1 ± 0.2	2%	86%	22.8 ± 0.3	3%	86%
	100	105.5 ± 2.3	4%	95%	115.5 ± 1.0	2%	85%
	200	213.3 ± 3.1	3%	93%	218.3 ± 3.7	3%	91%
PGE ₂	5	5.5 ± 0.1	4%	90%	4.9 ± 0.2	6%	97%
	20	18.0 ± 0.5	4%	90%	22.9 ± 0.9	7%	86%
	100	116.8 ± 4.0	7%	83%	113.8 ± 2.3	4%	86%
	200	206.5 ± 1.2	1%	97%	209.5 ± 4.6	4%	95%
PGF _{2α}	5	5.5 ± 0.1	2%	90%	5.2 ± 0.1	3%	97%
	20	19.2 ± 0.5	4%	96%	25.6 ± 0.4	3%	72%
	100	120.3 ± 1.7	3%	80%	123.8 ± 1.1	2%	76%
	200	233.3 ± 3.1	3%	83%	220.8 ± 0.9	1%	90%
TXB ₂	5	4.7 ± 0.2	7%	94%	5.0 ± 0.01	0.3%	100%
	20	16.8 ± 0.8	8%	84%	22.2 ± 1.3	10%	89%
	100	111.2 ± 4.1	7%	89%	118.0 ± 4.5	8%	82%
	200	219.0 ± 3.6	3%	91%	189.5 ± 6.9	7%	95%

Table 8 (part 3/4): Precision and accuracy of LM metabololipidomics. %RSD, relative standard deviation; %RSD = (SEM/mean) x 100; mean ± SEM; n = 4

Lipid Mediator Metabololipidomics *LC-MS-MS workup*



© Serhan et al.

Compound name	Nominal concentration (pg)	Intraday (repeatability)				Interday (reproducibility)			
		Measured (pg)		%RSD	Accuracy	Measured (pg)		%RSD	Accuracy
15-HETE	5	4.1 ± 0.1	2%	82%	4.2 ± 0.05	2%	84%		
	20	16.6 ± 0.2	2%	83%	21.8 ± 0.4	3%	91%		
	100	109.3 ± 2.6	5%	91%	110.8 ± 1.3	2%	89%		
	200	212.8 ± 2.3	2%	94%	219.0 ± 0.8	1%	91%		
12-HETE	5	4.5 ± 0.04	2%	89%	4.8 ± 0.04	2%	95%		
	20	17.7 ± 0.3	2%	89%	23.1 ± 0.6	4%	85%		
	100	113.8 ± 1.8	3%	86%	117.0 ± 1.1	2%	83%		
	200	216.5 ± 1.8	2%	92%	220.0 ± 1.7	2%	90%		
5-HETE	5	4.9 ± 0.05	2%	97%	5.0 ± 0.1	4%	100%		
	20	20.3 ± 0.3	3%	98%	25.7 ± 0.3	2%	72%		
	100	125.0 ± 3.0	5%	75%	127.8 ± 1.4	2%	72%		
	200	232.5 ± 1.4	1%	84%	239.5 ± 4.7	4%	80%		

Table 8 (part 4/4): Precision and accuracy of LM metabololipidomics. %RSD, relative standard deviation; %RSD = (SEM/mean) x 100; mean ± SEM; n = 4


Lipid Mediator Metabololipidomics LC-MS-MS workup



© Serhan et al.

	Day-0	Day-1	Day-2	Day-3	Day-4	Day-5	2-weeks	4-weeks
DHA bioactive metabolome		Δ change (%)						
RvD1	49.5±1.3							
with		-1.0	-0.2	+3.2	-3.2	-1.1		
without		+0.8	-2.0	+2.4	-13.2	-15.5	-33.4	-56.8
RvD2	13.4±0.8							
with		-1.3	+1.1	+6.4	+2.3	+1.1		
without		-9.2	+2.5	+2.2	-1.2	-19.2	-39.5	-59.0
RvD3	68.5±7.5							
with		-3.1	+8.0	+6.6	-4.4	+6.0		
without		-6.8	-8.5	-3.3	-11.3	-12.6	-18.9	-25.4
RvD5	24.9±2.2							
with		+6.5	-10.5	-18.0	-8.4	-16.8		
without		-16.4	-14.8	-15.1	-20.2	-17.0	-38.7	-44.0
RvD6	21.6±2.0							
with		+1.7	+2.2	-1.1	+0.4	-11.8		
without		+4.6	-9.8	-12.5	-17.2	-14.3	-19.5	-85.3
PD1	6.1±0.1							
with		-3.0	+5.5	-8.3	-11.9	-10.4		
without		+6.7	+1.7	-18.3	-26.8	-23.7	-33.0	-36.7
MaR1	9.0±0.8							
with		+2.3	+0.9	-8.0	-6.5	-17.8		
without		-6.5	-2.5	-2.6	-15.7	-24.7	-67.9	-73.1
EPA bioactive metabolome								
RvE1	13.1±0.9							
with		-3.3	-9.0	-13.9	-9.4	-13.1		
without		+7.6	-4.2	-22.0	-23.7	-24.6	-50.7	-57.1
RvE3	686.3±51.0							
with		-6.4	-7.6	-5.1	-7.2	-4.4		
without		-1.3	-5.2	-4.5	-3.5	-18.7	-40.9	-42.2
AA bioactive metabolome								
LXA ₄	11.2±1.0							
with		1.8	+0.5	-8.6	-8.3	-5.9		
without		4.4	-3.8	-5.2	-5.7	-3.8	-14.4	-39.9
LXB ₄	117.6±2.6							
with		-1.7	-3.6	-1.9	+5.2	-5.5		
without		-1.3	-5.8	-3.6	-5.8	-11.3	-17.7	-57.1
PGD ₂	55.1±1.9							
with		+9.8	-6.8	+3.3	+3.7	-15.5		
without		+7.6	-2.8	+5.8	+2.1	-16.9	-10.9	-26.7
PGE ₂	199.2±10.2							
with		+3.7	-1.2	+0.1	-11.8	-15.6		
without		-6.2	-12.2	-14.3	-20.8	-15.3	-11.9	-29.8
PGF _{2a}	143.3±2.5							
with		-0.4	-0.2	+1.2	+8.9	-13.6		
without		+8.0	+1.1	+1.6	+4.8	-9.0	-15.6	-34.1
TxB ₂	6067.0±128.4							
with		-5.2	-7.8	+7.4	-5.5	-12.5		
without		+9.2	-3.2	+7.3	-16.6	-19.2	-17.8	-34.2

Table 9: Human serum (250µl) was stored under -30°C with or without 1ml methanol. On the designated day (Day-0 to 4-weeks), LC-MS-MS based LM metabololipidomics were carried out. Results are expressed as pg/ml serum; mean ± SEM of d = 3 determinations and Δ change (%) from Day-0.



CETRI

Center for Experimental Therapeutics
and Reperfusion Injury

Multiple Reaction Monitoring Transitions

Lipid Mediator Metabololipidomics LC-MS-MS workup



© Serhan et al.


Compound	Q1	Q3	Compound	Q1	Q3
AA	303.3	259.1	15-HEPE	317.2	219.1
PGD ₂	351.3	233.2	18-HEPE	317.2	259.1
PGE ₂	351.3	175.1			
PGF _{2α}	353.3	193.1	DHA	327.3	283.1
TXB ₂	369.3	169.1	RvD1	375.2	215.2
LTB ₄	335.3	195.2	RvD2	375.2	141.2
LXA ₄	351.2	235.1	RvD3	375.2	147.1
LXB ₄	351.2	221.1	RvD5	359.2	199.1
5,15-diHETE	335.2	115.2	RvD6	359.2	159.1
5-HETE	319.2	115.1	MaR1	359.2	250.1
12-HETE	319.2	179.1	4,14-diHDHA	359.2	221.1
15-HETE	319.2	219.1	PD1	359.2	153.1
			4-HDHA	343.2	101.1
EPA	301.3	257.1	7-HDHA	343.2	141.1
LXA ₅	349.2	215.1	14-HDHA	343.2	205.1
LXB ₅	349.2	221.2	17-HDHA	343.2	245.1
5,15-diHEPE	333.2	235.2			
RvE1	349.2	195.1	d ₄ -PGE ₂	355.3	193.2
RvE2	333.3	253.2	d ₅ -LXA ₄	356.3	115.2
RvE3	333.3	201.1	d ₄ -LTB ₄	339.3	197.2
5-HEPE	317.2	115.1	d ₈ -5S-HETE	327.3	116.1
12-HEPE	317.2	179.1	d ₅ -RvD2	380.2	141.2

Table 1: MRM transition for LC-MS-MS based metabololipidomics of AA, EPA and DHA derived lipid mediators.

Compound	Q1	Q3	Compound	Q1	Q3
AT-LXA ₄	351.2	235.1	AT-RvD3	375.2	147.1
AT-LXB ₄	351.2	221.1	AT-RvD5	359.2	199.1
			AT-RvD6	359.2	159.1
AT-RvD1	375.2	215.2			
AT-RvD2	375.2	141.2	AT-PD1	359.2	153.1

Table 2: MRM transition for LC-MS-MS based metabololipidomics of AA and DHA derived aspirin-triggered lipid mediators.

Lipid Mediator Metabololipidomics *LC-MS-MS workup*



© Serhan et al.

Compound	Q1	Q3	Compound	Q1	Q3
20-OH-LTB ₄	351.2	195.1	Dihydro-RvD1	377.2	279.1
20-COOH-LTB ₄	365.2	195.1	7-oxo-RvD2	373.2	260.1
15-oxo-LXA ₄	349.2	233.1	16-oxo-RvD2	373.2	260.1
13,14-dihydro-LXA ₄	353.2	253.1	22-OH-RvD2	391.2	249.1
5-oxo-LXB ₄	349.2	219.1	Dihydro-RvD2	377.2	279.1
			14-oxo-MaR1	357.2	221.1
18-oxo-RvE1	347.2	213.1	7-oxo-MaR1	357.2	219.1
19-OH-RvE1	365.2	205.1	22-OH-MaR1	375.2	177.1
20-OH-RvE1	365.2	291.1	22-COOH-MaR1	389.2	177.1
20-COOH-RvE1	379.2	237.1	Dihydro-MaR1	361.2	223.1
Dihydro-RvE1	351.2	225.1	22-OH-PD1	375.2	153.1
			22-COOH-PD1	389.2	153.1
8-oxo-RvD1	373.2	231.1	Dihydro-PD1	361.2	263.1
17-oxo-RvD1	373.2	231.1			

Table 3: MRM transition for LC-MS-MS based metabololipidomics of further metabolites of AA, EPA and DHA derived lipid mediators.


Additional references:

Specific lipid mediator signatures of human phagocytes: microparticles stimulate macrophage efferocytosis and pro-resolving mediators. *Dalli J, Serhan CN. Blood.* 2012 120(15)


Metabolomics-lipidomics of eicosanoids and docosanoids generated by phagocytes. *Yang R, Chiang N, Oh SF, Serhan CN. Curr Protoc Immunol.* 2011 Chapter 14:Unit 14.26.

Prof. Charles N. SERHAN Laboratory
Brigham & Women's Hospital
Center for Experimental Therapeutics
and Reperfusion Injury
Harvard Medical School

Lipid Mediator Metabolipidomics ***LC-MS-MS workup***



© Serhan et al.



CNSerhanlab


CETRI

Center for Experimental Therapeutics
and Reperfusion Injury

Lipid Mediator Biosynthetic Pathways




Biosynthetic scheme for the Prostaglandins, Leukotrienes and Lipoxins




Key eicosanoids that play pivotal roles in initiating inflammation and its resolution.

○ Circle denotes human lipoxygenases, in mice 12/15 LOX is a single enzyme NOT found in human tissue.

Anti-Inflammatory and Pro-Resolving Lipid Mediators. *Charles N. Serhan, Stephanie Yacoubian and Rong Yang Annu Rev Pathol.* 2008 3: 279–312. [PMCID: PMC2739396](https://pubmed.ncbi.nlm.nih.gov/18711111/)




Biosynthetic scheme for the Leukotrienes



Leukotrienes and lipoxins: structures, biosynthesis, and biological effects. *Samuelsson B, Dahlén SE, Lindgren JA, Rouzer CA, Serhan CN Science. 1987 237(4819):1171-6 PMID:2820055*


Biosynthetic scheme for the E-series Resolvins



E-series resolvin biosynthesis and metabolic inactivation.

Resolvins and Protectins in Inflammation-Resolution *Charles N. Serhan, Nicos A. Petasis Chem Rev. 2011 12 111(10): 5922–5943. PMID: PMC3192290*

Biosynthetic scheme for the D-series Resolvins




Pathways and enzymes in the biosynthesis of D-series resolvins.

Resolvins and Protectins in Inflammation-Resolution *Charles N. Serhan, Nicos A. Petasis Chem Rev. 2011 12 111(10): 5922–5943. PMID: PMC3192290*



Biosynthetic scheme for the D-series Resolvins




Resolvin D3 and aspirin-triggered resolvin D3 are potent immunoresolvents *Dalli J, Winkler JW, Colas RA, Arnardottir H, Cheng CY, Chiang N, Petasis NA, Serhan CN Chem Biol. 2013 20(2):188-201. PMID: 23438748*

Resolvin D4 stereoassignment and its novel actions in host protection and bacterial clearance *Winkler JW, Orr SK, Dalli J, Cheng CY, Sanger JM, Chiang N, Petasis NA, Serhan CN Sci Rep. 2016 6:18972. PMID: 26743932*

Resolvin D3 multi-level proresolving actions are host protective during infection *Paul C. Norris, Hildur Arnardottir, Julia M. Sanger, David Fichtner, Gregory S. Keyes, Charles N. Serhan PLEFA. 2016 In Press, Accepted Manuscript*


Biosynthetic scheme for the aspirin triggered D-series Resolvins



Pathways and enzymes involved in the biosynthesis of aspirin-triggered D-series resolvins.

Resolvins and Protectins in Inflammation-Resolution *Charles N. Serhan, Nicos A. Petasis Chem Rev.* **2011** 12 111(10): 5922–5943. [PMCID: PMC3192290](https://pubmed.ncbi.nlm.nih.gov/2192290/)

Biosynthetic scheme for the Protectins




Resolvins and Protectins in Inflammation-Resolution *Charles N. Serhan, Nicos A. Petasis* *Chem Rev.* **2011** 12 111(10): 5922–5943. [PMCID: PMC3192290](https://pubmed.ncbi.nlm.nih.gov/23192290/)

Novel proresolving and tissue-regenerative resolvins and protectin sulfido-conjugated pathways *Dalli J, Ramon S, Norris PC, Colas RA, Serhan CN* *FASEB J.* **2015** 29(5):2120-36 [PMCID: 25713027](https://pubmed.ncbi.nlm.nih.gov/25713027/)

Identification and Actions of A Novel Third Maresin Conjugate in Tissue Regeneration: MCTR3. *Dalli J, Sanger JM, Rodriguez AR, Chiang N, Spur BW, Serhan CN.* **PLOS ONE** **2016**; in press.



Biosynthetic scheme for the Maresins



Maresins: novel macrophage mediators with potent antiinflammatory and proresolving actions. *Serhan CN, Yang R, Martinod K, Kasuga K, Pillai PS, Porter TF, Oh SF, Spite M. J Exp Med.* 2009 206(1):15-23. [PMCID: PMC2626672](https://pubmed.ncbi.nlm.nih.gov/19226672/)

Novel proresolving and tissue-regenerative resolvin and protectin sulfido-conjugated pathways *Dalli J, Ramon S, Norris PC, Colas RA, Serhan CN FASEB J.* 2015 29(5):2120-36 [PMCID: 25713027](https://pubmed.ncbi.nlm.nih.gov/25713027/)

The Protectin PCTR1 enhances resolution of infectious inflammation and is produced by human M2 macrophages. *Ramon S, Dalli J, Sanger JM, Winkler JW, Aursnes M, Tungen JE, Hansen TV, Serhan CN. Am J Pathol.* 2016; in press.